



## Mounting Suggestions: Honda VF750 Magna 1994-03 Tank Bib

### Mounting of Mustang Tank Bib

1. Remove seat and gas cap.
2. Align hole in tank bib with gas filler opening.
3. Secure front of tank bib with supplied Velcor® to the underside of the tank.
4. Remove rear tank bolt.
5. Align hole in bib with tank mount.
6. Secure with stock rear tank bolt.
7. Reinstall seat.



**Enjoy the custom look of your Mustang Tank Bib!**



These steps are guidelines to help you install your new Mustang product. If you need further assistance please call us at 800-243-1392, Monday through Friday 9:00 - 5:30 Eastern or send an e-mail to [questions@mustangseats.com](mailto:questions@mustangseats.com) or visit our website at <http://www.mustangseats.com>